

REUSE/Cleaning of Eye Protection

4/23/20

- When caring for COVID positive or suspected patients, Airborne/Contact/Eye Protection precautions are required.
- In effort to conserve use of PPE- Eye Protection, reuse and cleaning are measures we have implemented until supply chains are restored.
- Please follow the below steps to ensure proper cleaning of eye protection. Please be judicious in the reuse of these types of equipment.

Eye Protection:

- Eye Protection (goggles or face shield) can be cleaned using an A3 or Microkill wipe after each use



- Wipe down the entire apparatus
- Observe proper wet times as per manufacturer's recommendations
- **If the integrity of the Eye Protection is compromised – the viewing field is scratched making it difficult to see, it no longer fits tightly to your head or no longer remains in place when it is on, or any other concerns arise, please replace the Eye Protection immediately.**